

Incubator IN and IF with SingleDISPLAY
Incubator INplus and IFplus
with TwinDISPLAY
Natural convection or forced air circulation
AtmoCONTROL software

Model sizes: 30 / 55 / 75 / 110 / 160 / 260 / 450 / 750 +20 °C to +80 °C

**INCUBATOR I** Memmert incubators I are at home in the world of research, pharmaceutics and food technology. The heating and control system are especially optimised for low temperatures of up to +80 °C. To prevent temperature overshoots, temperature is increased within a very narrow control range and kept exactly at the setpoint value. As required, the models with natural convection or with forced air circulation are available.





## As little air circulation as possible in the incubator


Forced air convection may destroy the protective layer from moist air that is generated during incubation over the samples. This would lead to dehydration of the culture. In a Memmert incubator, the perfect combination of all-round surface heating and temperature control system ensures that incubation generally takes place without forced air circulation. Provided the chamber is fully loaded and forced air circulation is required, it can be precisely adjusted in 10 % steps from 0 to 100 %.

## Sterilisation

The chamber of the incubators INplus/IFplus, including all installations and sensors, can be sterilised at +160 °C in a 4-hour programme to guarantee optimum hygiene.

## Fresh air is preheated

Temperature deviations caused by fresh air can influence sample characteristics or prolong drying. In Memmert incubators, the fresh air is therefore fed through a pre-heating chamber and seamlessly introduced into the working chamber.



Air supply from outside

according to DIN 12880:2007-05, EN 61010-1 (IEC 61010-1), EN 61010-2-010







Interior: Stainless steel, material 1.4301 (ASTM 304) with

all-round deep-drawn ribs to integrate the large-area heating with ceramic-metal sheath

Textured stainless steel, rear zinc-plated steel, intuitively operated SingleDISPLAY or TwinDISPLAY (TFT colour display) with touchscreen; inner glass door, outside fully insulated stainless steel door (from size 450 two leaves)

Admixture of pre-heated fresh air by electronically adjustable airflap Fresh air:

Mains cable with plug (German type) Connection:

Installation: 4 feet; sizes 450/750 mounted on lockable castors


Interfaces:

Housing:





USB: only TwinDISPLAY



Model sizes/Descri	ption		30	55	75	110	160	260	450	750		
Stainless steel interior	Volume	approx. l	32	53	74	108	161	256	449	749		
	Width	(A) mm	400		560		640	10	40			
	Height	(B) mm	320	400	560	480	720	800	720	1200		
	Depth (less 39 mm for fan)	(C) mm	250	250 330		400		500	60	00		
	Max. number of grids/shelves	number	3	4	6	5	8	9	8	14		
	Max. loading per grid/shelf	kg	20 30									
	Max. loading of chamber	kg	60					300				
	Max. loading per slide-in drip tray	kg	1,5				3	4	8			
	Max. loading per bottom drip tray	kg	1,5				3	4	8			
Textured stainless steel exterior	Width	(D) mm		585		745		824	1224			
	Height (size 450, 750 with castors)	(E) mm	704	784	944	864	1104	1183	1247	1720		
	Depth (without door handle, depth of handle +56 mm)	(F) mm	434	5	14	58	34	684	784			
Standard equipment	Stainless steel grids, electropolished	number	1 2									
	Standard works calibration certificate (measuring point chamber center)	°C	+37									
Temperature	Working temperature range	°C	at least 5 (IN/INplus) 10 (IF/IFplus) above ambient temperature up to +80									
	Setting temperature range	°C	+20 to +80									
	Setting accuracy	°C	0.1									
Further data	Electrical load at 230 V, 50/60 Hz	approx. W	1600	1000	1250	1400	1600	1700	1800	2000		
	Electrical load at 115 V, 50/60 Hz	approx. W	800	800 900 1500						1800		
Packing data	Net weight	approx. kg	48	57	66	76	96	110	161	217		
	Gross weight (packed in carton)	approx. kg	64	76	85	101	122	161	227	288		
	Width	approx. mm	660	730		830		930	1330			
	Height	approx. mm	890	950	1130	1050	1300	1380	1440	1910		
	Depth	approx. mm	650	6	70	80	00	930	1050			
Order No. Incubators		IN30	IN55	IN75	IN110	IN160	IN260	IN450	IN750			
I = Incubator		IN30plus	IN55plus	IN75plus	IN110plus	IN160plus	IN260plus	IN450plus	IN750plus			
N = Natural convection		IF30	IF55	IF75	IF110	IF160	IF260	IF450	IF750			
F = Forced convection plus = Model with TwinDISPLAY		IF30plus	IF55plus	IF75plus	IF110plus	IF160plus	IF260plus	IF450plus	IF750plus			

Accessories	30	55	75	110	160	260	450	750	
Stainless steel grid, electropolished		E20164		E20165		E28891	E20182		
Reinforced stainless steel grid, electropolished, max. loading 60 kg; from size 450 with guide bars and fixing screws (only in connection with option K1). Please consider max. loading of chamber		-		E29767		E29766	B32190		
Perforated stainless steel shelf	B29727 B03916		B00325		B29725	B00328			
Reinforced stainless steel shelf, max. loading 60 kg; with guide bars and fixing screws (requires option K1). Please consider max. loading of chamber		-						B32191	
Stainless steel slide-in drip tray, 15 mm rim (may affect the temperature distribution, not in connection with option K1)	E02070	E02070 E02072		E02073		E29726	26 E02075		
Stainless steel slide-in drip tray, 15 mm rim, with guide bars and fixing screws (may affect the temperature distribution, only in connection with option K1)	-							B32763	
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution, not in connection with option K1)	B04356	66 B04358		B04359		B29722 B0436		362	
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution, only in connection with option K1)	-					B34055			
Wall bracket for wall mounting	B29755	B29756	B29757	B29758	B29759		-		
Guarantee extension by 1 year			GA2Q5						
USB-Ethernet adapter	E06192								
Ethernet connection cable 5 m for computer interface				E06189					

Accessories	30	55	75	110	160	260	450	750		
USB User-ID stick (with User-ID licence): Oven-linked authorisation licence (User-ID-programme) on Memory-stick, prevents undesired manipulation by unauthorised third parties. When reordering please specify serial number (only for units with TwinDISPLAY)		B33170								
USB stick with documentation software AtmoCONTROL and operation manual for products with SingleDISPLAY (the standard equipment of appliances with TwinDISPLAY includes one USB stick with AtmoCONTROL). When reordering please specify serial number	B33172									
et of height adjustable feet (4 pcs)		B29768					-			
tacking set (4 pcs) for stacking of appliances of same size		B29744					-			
Plug-in tube extension (outer diam. 60.3 mm, inner 57 mm), straight, for exhaust air ducting (if necessary for connection by hose)	B29718									
Plug-in tube extension (outer diam. 60.3 mm, inner 57 mm), angled, for exhaust air ducting (if necessary for connection by hose)	B29719									
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), with air slots	B29728	B29730	B29732	B29734	B29736	B29738	B29740	B29742		
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), without air slots	B29729	B29731	B29733	B29735	B29737	B29739	B29741	B29743		
Subframe, adjustable in height (size 30 to 75: height 600 mm, size 110 to 450: height 500 mm)	B29745	B29747		B29749		B29751	B29753	-		
Subframe, on castors (size 30 to 75: height 660 mm, size 110 to 160: height 560 mm)	B29746	B29	B29748 B2		9750		-			
subframe, adjustable in height, height 130 mm, for example for units with fresh air filter		B33657 B33659		B33661		B33664 -		-		
Software conforming to FDA AtmoCONTROL. Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Base licence for the control of one unit (only for units with TwinDISPLAY). Respective IQ/OQ documents available in German and English language (without surcharge)	FDAQ1									
Integration of one additional unit (up to max. 31 units) into an already existent FDA-software licence (only for units with TwinDISPLAY)	FDAQ2									
DAkkS calibration for one free-selectable temperature value according to method C (DKD-R 5-7)		E39696								
DAkkS calibration for further temperature values according to method C (DKD-R 5-7)		E39697								
IQ document with device-specific works test data, OQ/PQ check list as support for validation by customer		D00124								
IQ/OQ document with device-specific works test data for one free-selectable temperature value, incl. temperature distribution survey at Memmert for 9 measuring points (size 30), 27 measuring points (sizes 55 - 1060) to DIN 12880:2007-05. PQ check list as support for validation by customer. 305 € for further temperature values	D00125 D00127									
On-site IQ/OQ for a freely selectable temperature value, including temperature distribution survey for 9 measuring points (size 30), 27 measuring points (sizes 55 - 750) to DIN 12880: 2007-05 (excluding travel costs, not subject to discount, GER, AT, FR only)	DLQ100			2100						
Extension of DLQ100 by an additional freely selectable temperature value (not subject to discount)	DLQ100A									
Individual on-site Performance Qualification (PQ)	DLQ200									
Maintenance UIS - carrying out and documentation according to Memmert maintenance plan (excluding travel costs, not subject to discount, GER, AT, FR only)	S00311									
Maintenance contract UIS - carrying out and documentation according to Memmert maintenance plan, minimum duration 3 years (excluding travel costs, not subject to discount, GER, AT, FR only)	S00311J									
Calibration of one freely selectable temperature value (excluding travel costs, not subject to discount, GER, AT, FR only)	S00205									
Calibration of an additional temperature value (not subject to discount)  S00215										